## IMMOBILIZATION OF GLUCOAMYLASE FOR THE HYDROLYSIS OF STARCH AND DEXTRINIZED STARCH

Department of Biochemistry, Faculty of Medicine,
University of Jaffna

The isoelectric pH of glucoamylase from Aspergillus niger was between 3.5 and 4.0. In the studies for the selection of carrier for immobilization of glucoamylase cation exchanger Amberlite IR-120 and anion exchanger Amberlite IRA-904 were used. For Amberlite 1R-120, immobilization yields of protein and glucoamylase units were 10 and 11% respectively at pH 3.3 and at pH 2.5 were 14 6 and 15.1% respectively. Retentivities of Amberlite IR-120 for glucoamylase activity at 55°C and at pH 3.3 and 2.5 after 5 successive batch reactions (each 10 min duration) were 1.6 (+/-0.47) and 2.7 (+/-0.47) respectively. Cation exchanger Amberlite IR-120 was not suitable for immobilization of this enzyme as it showed poor retaining capacity. Hence the experiments were carried out with anion exchanger Amberlite IRA-904. The capacity was 42mg protein (493 AMG units) per gram (wet Wt) of resin at pH 4.5, giving an immobilization yield of 64% protein and 67.5% AMG units. The activity yield (a measure of functional fraction of expected immbolized enzyme activity) was 2.2% which indicates poor expression of immobilized glucoamylase activity. The retentivity of the enzyme activity was 94.4% (+/ - 2.3). Relative activities of soluble glucoamylase with 2% (W/V) starch, dextrinized starch (DE 36) and maltose were 100, 52.2 and 37.5% respectively. The activity yields of Amberlite IRA-904 immobilized glucoamylase were 22.3 and 28.0% with dextrinized starch and maltose respectively. Effectiveness factor, which is a measure of diffusional restriction for substrate molecule was 7.3, 1.6 and 1.1 for starch, dextrinized starch and maltose respectively. Thus glucoamylase is suitable for the physical immobilization on Amberlite IRA-904 and for continuous saccharification. It is desirable to use dextrinized starch than soluble starch and maltose with the immobilized glucoamylase.